

ABGENEX Pvt. Ltd., E-5, Infocity, KIIT Post Office, Tel: +91-674-2720712, +91-9437550560 Email: info@abgenex.com

Bhubaneswar, Odisha - 751024, INDIA

## 20-1040: Polyclonal antibody to Caspase-3 (active/cleaved)

Clonality: Polyclonal Application: IP,IHC,WB

**Reactivity:** Dog,Rat,Mouse,Human

 Gene :
 CASP3

 Gene ID :
 836

 Uniprot ID :
 P42574

 Format :
 Sera

Alternative Name: Apopain, Cysteine protease CPP32, Protein Yama, SREBP cleavage activity 1, CPP32

**Isotype:** Rabbit IgG

Immunogen Information: A recombinant catalytically active human caspase-3 protein was used as immunogen for this

antibody

## **Description**

Apoptosis, or programmed cell death, is a common property of all multicellular organisms. The current dogma of apoptosis suggests that the components of the core cell-death machinery are integral to cells and widely conserved across species. Caspases, a family of cysteinyl aspartate-specific proteases, are integral components of the cell death machinery (reviewed in Siegal, 2006; and Lavrik et al, 2005). They play a central role in the initiation and execution of apoptotic cell death and in inflammation. Caspases are typically divided into 3 major groups, depending on the structure of their prodomain and their function. Group 1: inflammatory caspases (caspases 1, 4, 5, 11, 12, 14). Group II: initiator of apoptosis caspases (caspases 2, 8, 9). Group II: effector caspases (caspases 3, 6, 7). Caspases are synthesized as zymogens (inactive pro enzyme precursors which require a biochemical change to become active enzymes) with an N-terminal prodomain of variable length followed by a large subunit (p20) and a small subunit (p10). Caspases are activated through proteolytic cleavage at specific asparagine residues that are located within the prodomain, the p10, and p20 subunits. Activation results in the generation of mature active caspases that consist of the heterotetramer p202-p102. Active caspases mediate cell death and inflammation through cleavage of particular cellular substrates that are involved in these processes. The Caspase-3 polyclonal antisera recognizes the proform of caspase-3 (approx. 32) kDa), and the large (approx. 14-21 kDa) and small (approx. 10 kDa) subunits of active/cleaved Caspase-3.

## **Product Info**

**Amount :** 50 μl **Content :** 50 μl sera

Storage condition:

Storage condition:

Store the antibody at 4°C, stable for 6 months. For long-term storage, store at -20°C. Avoid

repeated freeze and thaw cycles.

## **Application Note**

WB: 1:1000-1:2000, IHC (paraffin): 1:1000-1:5000, IHC (frozen): Users should optimize, IP: 1:50-1:200



ABGENEX Pvt. Ltd., E-5, Infocity, KIIT Post Office, Tel: +91-674-2720712, +91-9437550560 Email: info@abgenex.com

Bhubaneswar, Odisha - 751024, INDIA




Fig:1 Formalin-fixed, paraffin-embedded section of dog ischemic brain stained for Active/Cleaved Caspase-3 expression using 20-1040 at 1:2000. Staining is seen in the nuclei of dying neurons (black arrow) but not in the morphologically normal nuclei (open arrows). Caspase-3 expression in the nucleus is considered to be a marker of active/caspase-3 expression and apoptosis. Hematoxylin-eosin counterstain.




Fig:2 Formalin-fixed, paraffin-embedded section of irradiated mouse spleen stained for Active/Cleaved Caspase-3 expression using 20-1040 at 1:2000. Staining is seen in the nuclei of a subset of the cell population. Caspase-3 expression in the nucleus is considered to be a marker of active/caspase-3 expression and apoptosis. Hematoxylineosin counterstain.