

ABGENEX Pvt. Ltd., E-5, Infocity, KIIT Post Office, Tel: +91-674-2720712, +91-9437550560 Email: info@abgenex.com

Bhubaneswar, Odisha - 751024, INDIA

## 32-2203: CA11 Recombinant Protein

Alternative Name: Carbonic Anhydrase XI, Carbonic Anhydrase-Related Protein 2, Carbonic Anhydrase-Related Protein 11, CARP-2, CA-RP II, CARP XI, CARPX1, CA-XI.

## **Description**

Source: Escherichia Coli. CA11 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 326 amino acids (24-328) and having a molecular mass of 36.3kDa.CA11 is fused to a 21 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques. Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes which catalyze the reversible hydration of carbon dioxide. These metalloenzymes participate in various biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. The metalloenzymes exhibit extensive diversity in tissue distribution and in their subcellular localization. Carbonic Anhydrase XI (CA11) is probably a secreted protein, nevertheless, drastic changes at active site residues completely conserved in CA isozymes with catalytic activity, make it unlikely that CA11 has carbonic anhydrase activity. CA11 shares properties in common with 2 other acatalytic CA isoforms, CA VIII and CA X. CA11 is amply expressed in the brain, and may have a general role in the central nervous system.

## **Product Info**

Amount: 20 µg

**Purification:** Greater than 85% as determined by SDS-PAGE.

Content: The CA11 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.4M Urea and 10% glycerol.

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of

Storage condition: time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid

multiple freeze-thaw cycles.

Amino Acid: MGSSHHHHHH SSGLVPRGSH MHIGPAPDPE DWWSYKDNLQ GNFVPGPPFW GLVNAAWSLC

AVGKRQSPVD VELKRVLYDP FLPPLRLSTG GEKLRGTLYN TGRHVSFLPA PRPVVNVSGG PLLYSHRLSE LRLLFGARDG AGSEHQINHQ GFSAEVQLIH FNQELYGNFS AASRGPNGLA ILSLFVNVAS TSNPFLSRLL NRDTITRISY KNDAYFLQDL SLELLFPESF GFITYQGSLS TPPCSETVTW ILIDRALNIT SLQMHSLRLL SQNPPSQIFQ SLSGNSRPLQ PLAHRALRGN

RDPRHPERRC RGPNYRLHVD GVPHGR

