

ABGENEX Pvt. Ltd., E-5, Infocity, KIIT Post Office, Tel: +91-674-2720712, +91-9437550560 Email: info@abgenex.com Bhubaneswar, Odisha - 751024, INDIA

32-5657: Recombinant Hemagglutinin-Influenza A Virus H3N2 Perth 16/09

Description

Source: Baculovirus H3N2 produced in Hi-5 cell of Baculovirus is a single polypeptide chain containing 339 amino acids (17-345) and having a molecular mass of 37.8kDa.H3N2 is fused to a 6 amino acid His-tag at C-terminus & purified by proprietary chromatographic techniques. H3N2 is a subtype of the influenza A virus. Its name derives from the forms of the two kinds of proteins on the surface of its coat, hemagglutinin (H) and neuraminidase (N). H3N2 exchanges genes for internal proteins with other influenza subtypes. H3N2 has tended to dominate in prevalence over H1N1, H1N2, and influenza B. H3N2 strain descended from H2N2 by antigenic shift, in which genes from multiple subtypes re-assorted to form a new virus. Both the H2N2 and H3N2 strains contained genes from avian influenza viruses.H3N2 viruses are able to infect mammals and birds. In pigs, humans, and birds, the virus has mutated into many strains. Hemagglutinin(HA) binds to sialic acid-containing receptors on the cell surface, generating the attachment of the virus particle to the cell. HA has a vital part in the determination of host range restriction and virulence and is in charge of the diffusion of the virus into the cell cytoplasm by facilitating the fusion of the membrane of the endocytosed virus particle with the endosomal membrane.

Product Info

Amount: 5 µg

Purification: Greater than 90% as determined by SDS-PAGE.

Content: The H3N2 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0) and 10% glycerol.

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of Storage condition:

time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid

multiple freeze-thaw cycles.

Amino Acid: ADPMQKLPGN DNSTATLCLG HHAVPNGTIV KTITNDQIEV TNATELVQSS STGEICDSPH

> QILDGKNCTL IDALLGDPQC DGFQNKKWDL FVERSKAYSN CYPYDVPDYA SLRSLVASSG TLEFNNESFN WTGVTQNGTS SACIRRSKNS FFSRLNWLTH LNFKYPALNV TMPNNEQFDK LYIWGVHHPG TDKDQIFLYA QASGRITVST KRSQQTVSPN IGSRPRVRNI PSRISIYWTI VKPGDILLIN STGNLIAPRG YFKIRSGKSS IMRSDAPIGK CNSECITPNG SIPNDKPFQN

VNRITYGACP RYVKQNTLKL ATGMRNVPEK QTRHHHHHH.

